

Technical Specifications

Measure system	White LED and Color Sensor
Throughput	140 tests/hour(fast mode),50 tests/hour(normal mode)
Memory	1000 Recent Test Results
Printer	Wireless Thermal Printer(optional)
Interface	Mini USB
Keypad	Touch Keypad
Display	LCD Screen
Voltage	DC 3.7V
Adapter Input	AC 100-240V,50/60Hz, 2W
Continuous test time	>6 Hours
Power	10 VA
Net weight	160g
Dimension(mm)	110(L)*68(W)*27(H)



A hand is holding a white handheld urine analyzer. The device has a color LCD screen displaying test results for various parameters: NO. 001, TYPE 11, LEU -, BLD +1, NIT -, KET -, UBG -, BIL +2, PRO +2, GLU +3, pH 5.0, VC -, SG 1.020. The screen also shows a date field. Below the screen is a touch keypad with a power button, a back arrow, and four directional arrows. The device is labeled 'Urine Analyzer' on the side.

Ui

Handheld Urine Analyzer

New Favorite in IVD Doctors' Good Mobile Assistant

EMFSUN CHENGDU EMFSUN MEDICAL TECHNOLOGY CO., LTD.

Add: No.333,Hezuo Road,West Hi-tech Zone,Chengdu,Sichuan,China.

Domestic Sales Tel: +86-028-6758-5399

International Sales Tel: +86-028-6558-9558

Customers Service Tel: +86-400-605-9899

Website: www.empsun.com

All rights are reserved by Chengdu Empsun Medical Technology Co.,LTD.like changing product item number,specification,appearance,etc. Pls be subject to the latest leaflet.EN .UI.201808

DISTRIBUTOR



EN.UI.201808



EMFSUN

Global Expert in Rapid Diagnosis
Good Sale in Domestic and Foreign Market



Product Advantage



Patents Protected

To be made in national High-Tech company with several patents.



Bluetooth Technology

To connect by BT4.0 with low power dissipation, long distance transmit and no need matching password.



1000 Results Storage Capacity

To check patients' historical data easily by super storage capacity.



Easy and Noninvasive

Professional home use IVD device. Easy and noninvasive operation.



Only 160g. Fashion and Mini

The first smallest urine analyzer in the world. Appearance is fashion and delicate.



Wireless Thermal Printing

2.4Ghz wireless thermal printer is optional. It can print international unit report in Chinese or English.

User type



Nephrotic Patient

Can monitor the state of illness without going to hospital often.



Middle and Old Age People

Daily test to grasp the best time of cure.



Diabetic

Monitoring renal function to alert patient about complication.



People Under High Pressure

Testing immunity to adjust work status.



Expectant Mothers

Daily test to prevent the disease during pregnancy.



Infants

Routine test to detect the state of growth.

Test Item

Leukocytes(LEU)

Urobilinogen(UBG)

PH

Specific-gravity(SG)

Bilirubin(BIL)

Ketone(KET)

Nitrite(NIT)

Protein(PRO)

Blood(BLD)

Vitamin C(VC)

Glucose(GLU)

- Fashion and Mini Urine Analyzer
- Wireless Thermal Printer(optional)
- Built-in High-capacity Li-battery
- Precise and Efficient



Wireless printer